

MS-V 1 CCD ENDO-CAMERA

Order-No. 670-15000



Higher Contrast

The Endo-Camera *MS-V* can provide an enhanced resolution. Even fine tissue structures can be visualized with its detailed transmission.

Great Light Sensitivity

Under difficult light conditions the Endo-Camera *MS-V* is showing a well illuminated picture. For small diameter telescope's the **GAIN-Key** should be used.

New Camera Head-Design

The Endo-Camera *MS-V* is specially designed for endoscopic procedures. The angled camera-head gives a perfect sense of up and down to the user. The new camera is much lighter and smaller. The head button is individual selectable.

White Balance Memory Function

The white balance setting is stored in the camera's micro-processor. Thus, the correct adjustment is kept even if the Endo-Camera *MS-V* is switched off.

USB II out

With **USB II**, you can directly connect the *MS-V* to PC. This allows to see and store the video on PC.




Head Button

Individual selectable head button: "W/B", "Print" and "Gain" (switch on the backside).



Moiré Function

Allows to use flexible Fibrescopes.



Freeze

Screen Splitting

Freeze and Screen Splitting

MS-V allows – controlled by a footswitch – **freeze** and **screen splitting** on the monitor for detailed comparison.