

Dentistry, Dermatology



## *KLS Martin MD 62*

THE COMPACT ELECTROSURGICAL UNIT  
BUILT TO HIGH SAFETY STANDARDS

**KLS** martin  
GROUP

HF surgery in perfection

## KLS Martin MD 62 Electrosurgical Unit

*The KLS Martin MD 62 is an electrosurgical unit for use in dentistry.*

*Additional fields of application include dermatology and minor cosmetic surgery.*



### KLS Martin MD 62

HF surgery is considered a safe and proven method in dentistry. With its special features and accessories and matching output power, KLS Martin's MD 62 electrosurgical unit meets all requirements in this field.

Thanks to its small dimensions, the unit is very easy to handle, providing maximum flexibility in use.

The KLS Martin MD 62 has the following main features:

- Two different types of current, one for cutting and one for coagulation.
- To account for the special conditions and requirements associated with dental, oral and maxillomandibular applications, the high-frequency current can be activated in two different ways. In hand switch mode, the HF current is activated via a finger-operated switch located on the handle. In foot switch mode, a handle without a pushbutton is used and the HF current is activated via a foot switch.
- An optional foot switch is available for use with a handle without pushbutton.
- The powerful electrosurgical unit is supplemented by a wide range of high-quality electrodes and accessories.



### Safe and Easy Operation

All functions available are clearly marked by symbols on the front panel. The unit features an analog rotary switch for safe and easy power setting.

The two types of current can be activated – at the user's option – either via a foot switch or a handle featuring a pushbutton.

All connector jacks are marked unambiguously for correct use.

### Safety in Use

The KLS Martin MD 62 has been designed to meet the high demands of surgical users.

The unit bears the CE-mark in compliance with Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Due to the monopolar technique, no neutral electrode is needed. This spares the patient the inconvenience caused by NE attachment.

### Performance Characteristics

The user has two types of current available: one for cutting and one for coagulation. These currents have been matched to the tissue-specific impedances involved and are fully compatible with the monopolar mode of operation.

This powerful electrocautery device provides the high technical quality, reliability and safety medical specialists require in daily practice.

### Hygiene Made Easy

The housing is made of high-quality, coated sheet steel. This reliably prevents oxidation caused by disinfectants.

As neither cooling slots nor fans are required for safe operation, no liquids such as cleaning or disinfecting solutions can penetrate into the unit.



### Wide Range of High-Quality Accessories

Handles specially optimized for use in the oral cavity plus a comprehensive range of active electrodes guarantee the unit's usefulness and reliability in daily practice.

Selected accessories have been combined into a practical set.

### Application examples from dentistry:

- Gingivoplasty/gingivectomy
- Dentitis difficilis
- Hemostasis by using HF current (electrocoagulation)
- Biopsy
- Tumor resection/coagulation
- Root canal sterilization
- Mucosal surgery (removal of abscesses, overgrowth, fibromas)
- Lengthening of tooth crowns or abutment teeth that are too short for bridge support
- Surgical treatment of abnormal labial or lingual frenula
- Desensitization of bare dental necks

## KLS Martin Subsidiaries

### **KLS Martin France SARL**

3, Rue Gambetta  
68000 Colmar  
France  
Tel. +33 3 89 21 66 01  
Fax +33 3 89 23 65 14  
axel.siegelin@klsmartin.com

### **Martin Nederland/Marned B.V.**

Visserstraat 9a  
1271 VE Huizen  
P.O. Box 278  
1270 AG Huizen  
The Netherlands  
Tel. +31 35 523 45 38  
Fax +31 35 523 53 48  
info@martinnederland.nl

### **KLS Martin L.P.**

11239-1 St. John's Industrial  
Parkway South  
Jacksonville, FL 32246  
USA  
Office phone +1 904 641 7746  
Office fax +1 904 641 7378  
Toll free +1 800 625 1557  
sburke@klsmartin.com  
www.klsmartinusa.com

### **Martin Italia S.r.l.**

Via Paracelso, 18  
Palazzo Andromeda, Ingresso 2  
20041 Agrate Brianza (Mi)  
Italy  
Tel. +39 039 605 67 31  
Fax +39 039 605 67 42  
info@martinitalia.it

### **Nippon Martin K.K.**

6 Fl., 4-6-16 Hiranomachi  
541-0046 Chuo-ku,  
Osaka 541-0046  
Japan  
Tel. +81 6 62 28 90 75  
Fax +81 6 62 28 90 76  
info@n-martin.com

### **Orthosurgical Implants Inc.**

12244 SW 130 Street  
Miami, FL 33186  
USA  
Office phone +1 305 969 4545  
Office fax +1 305 969 4546  
Toll free +1 888 305 4559  
sales@orthosurgical.com  
www.orthosurgical.com

### **Gebrüder Martin GmbH & Co. KG**

A company of the KLS Martin Group  
Ludwigstaler Str. 132 · D-78532 Tuttlingen  
Postfach 60 · D-78501 Tuttlingen  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com