

## ELICA II: VIDEO-CAMERA COMBINED WITH DUO LIGHT-SOURCE

Order-No. 670-51500



### Higher Contrast

The Endo-Camera *ELICA II* can provide an enhanced resolution. Even fine tissue structures can be visualized with its detailed transmission.

### Great Light Sensitivity

Under difficult light conditions the Endo-Camera *ELICA II* is showing a well illuminated picture. For small diameter telescope's digital filters allow to enhance the sharpness.

### New Camera-Head-Design

The Endo-Camera *ELICA II* is specially designed for endoscopic procedures. The angled camera-head gives a perfect sense of up and down to the user.

### White Balance Memory Function

The white balance setting is stored in the camera's micro-processor. Thus, the correct adjustment is kept even if the Endo-Camera *ELICA II* is switched off.

### USB II out (option)

With USB II, you can directly connect *ELICA II* to a PC. This allows to see and store the video on PC.



### Head Button

Individual selectable head button:  
 ■ "W/B": White balance  
 ■ "Print": Remote to recorder device  
 ■ "Sharp": Sharpness (switch on the backside).



### Moiré Function

Allows to use flexible Fibrescopes.



### Brightness

Light-intensity regulation.



Freeze

Screen Splitting

### Freeze and Screen Splitting

*ELICA II* allows freeze and screen splitting on the monitor for detailed comparison.